




ACTIVITY@WHEEL

PHYSICAL MONITORING



HOUSING



TECNIPLAST

innovation through passionSM



Activity Wheel Physical Monitoring

1/Easy monitoring, easy maintenance

Designed for easy monitoring of rodents physical activity levels. The cage is easily demountable and each single part can be washed/decontaminated.

The cage features a high level of transparency which allows for good visibility and inspection of the animals.

2/Flexibility: adaptable for rats and mice

This rodent Activity Wheel provides an easy, convenient method for measuring lab rodents' physical activity in response to chemical or environmental stimuli.


It is especially usefull for researches involving circadian rhythms or pharmaceutical testing.

Regular activity keeps muscles coordination of movements in optimal conditions. It increases and manteins physical fitness. It helps the control of body weight.

This product is available both for mice and rats.

Easy monitoring, easy maintenance

1



Flexibility: adaptable for rats and mice

2



Options

- Individually ventilated Activity Wheel cage.

TECHNICAL SHEET

2154F0105 for Rat/hamster

Overall Dimensions (H x W x D) - mm: 480 x 315 x 470

1284L0105 for Mice

Overall Dimensions (H x W x D) - mm: 364 x 258 x 350

NOTE:

- Complete Activity Wheel System for mice also available for bigger cages (1291H and 1290D).
- Complete Activity Wheel System for rat also available for bigger cages (1500U).



TECNIPLAST

HEADQUARTERS

ITALY | Tel. +39 0332 809711 • www.tecniplast.it • E-mail: tecnicom@tecniplast.it

SISTER COMPANIES

AUSTRALIA/NEW ZEALAND | Tel. + 61 2 8845 6500 • www.tecniplast.it • E-mail: info@tecniplast.com.au

FRANCE | Tel. 04 72 52 94 41 • www.tecniplast.fr • E-mail: info@tecniplast.fr

GERMANY | Tel. 08805 921320 • www.tecniplast.de • E-mail: info@tecniplast.de

UNITED KINGDOM | Tel. 0845 0504556 • www.tecniplast.it • E-mail: info@tecniplastuk.com

JAPAN | Tel. + 03 5770 5375 • www.tecniplast.it • E-mail: info@tecniplastjapan.co.jp

USA/CANADA | Toll Free: 877.669.2243 • www.tecniplastusa.com • E-mail: info@tecniplastusa.com

To find your local distributor, please visit www.tecniplast.it



ISO-9001 - Cert. n° 0875



ISO-14001 - Cert. n° 0256A/0

